

DxC 700 AU & AU680 Contamination Prevention Parameters							
No.	PRECEDING TEST NAME	FOLLOWING TEST NAME	REAGENT PROBE CLEANER	WASH COUNT	EFFECTIVE OF WATER CLEANING	PREVENT USE	
						Mixer	Cuvette
1.	HCY B08176	Lip 030	CLN 100%*	5	No	Yes	No
2.	UIBC 205	HCY B08176	Water	1	Yes	Yes	No
3.	HCY B08176	UIBC 205	Water	1	Yes	Yes	No
4.	HCY B08176	LDH 026	Water	1	Yes	Yes	Yes
5.	HCY B08176	LDH 027	Water	1	Yes	Yes	Yes
6.	HCY B08176	LDH 028	Water	1	Yes	Yes	Yes

* CLN OSR0001/ODR2000, replenish cleaner on a daily basis or as requested.